



connect - share - transform

TMAX Kinesiology Tape (Pink)

Professional Muscle & Joint Support for Physiotherapy, Sports & Rehabilitation



TMAX Kinesiology Tape is designed to provide targeted support for muscles and joints without restricting natural movement. Suitable for athletes, physiotherapists and rehabilitation, it features breathable cotton fabric, waterproof acrylic adhesive and high elasticity for all-day comfort.

Key Features

- High stretch cotton fabric
- Skin-friendly waterproof adhesive
- Supports muscles & joints
- Allows full range of motion
- Suitable for sports, gym & physiotherapy
- Latex free

Specification	Details
Brand	TMAX
Product	Kinesiology Tape
Width	2 in (5 cm)
Color	Pink
Material	Cotton Fabric + Acrylic Adhesive
Waterproof	Yes
Breathable	Yes
Latex Free	Yes
Application	Muscle Support & Joint Stabilization
Suitable For	Physiotherapy, Sports, Gym, Running
Body Areas	Shoulder, Neck, Back, Arms, Elbows, Wrists, Knees, Ankles, Calves

Frequently Asked Questions

What is it used for?

Supports muscles, joints and rehabilitation.

Is it waterproof?

Yes.

Can I wear it during exercise?

Yes.

Which body parts?

Shoulders, back, knees, ankles and more.

Does it restrict movement?

No, it stretches naturally.

Sensitive skin?

Latex-free; patch test recommended.

Wear duration?

Typically 3–5 days.

Can it be cut?

Yes, into I, Y, X or fan strips.

Reusable?

No.

Who can use it?

Athletes, physiotherapists and fitness users.

Contact Us

MediKart Healthcare Systems Pvt. Ltd.

■ +91 7669011322

✉ info@medikart.co.in

■ www.medikart.co.in