

PORTABLE PHLEGM SUCTION UNIT

Product Code : MFPSM02

Oxy-med®
O₂ for healthy living



Specifications :

Input Power	AC 230V, 50Hz
Capacity	90VA
Negative Pressure Limit	> 0.08MPa
Pressure adjustment range	0.02MPa~Pressure Value Limit
Pumping Rate	> 15 L/min
Fuse	AC250V, F2.0 AL, ϕ 5 X 20
Reservoir bottle	1000ml, One Piece
Noise	< 60dB(A)
Vacuum gauge	0~0. IMP, 2.5 grade)
Equipment Type	High vacuum/low flow
Working mode	intermittent operation only 30 minutes ON, 30 minute OFF
Electrical requirements	is a class II, type BF

Note: It is not suitable for use in environment with flammable and explosive gases; Electromagnetic interference source: the minimum separation distance of 5s 1.1m when it is close to power lines and electrical equipment with the voltage less than 500V