

DreamStation Auto BiPAP features

- Advanced event detection*
- · Therapy modes: CPAP, fixed bi-level, auto bi-level
- · Adjustable PS min and PS max
- · NEW Smart Ramp (0-45 min)
- Flex Comfort: C-Flex, C-Flex+, Bi-Flex
- · **NEW** Performance Check
- Encore patient management system: compliance and data reports, waveforms reports
- **NEW -** Multiple remote monitoring solutions including integrated Bluetooth® technology and optional cellular or wifi modems
- · On-board memory
- · Oximetry module (optional)

Humidification

- $\boldsymbol{\cdot}$ Active humidification controlled by the room environment
- \cdot Pre-heat humidifier capable
- · Heated tube (optional)

DreamStation Auto BiPAP benefits


- Daily Progress Feedback displays the length of sleep and progress, designed to motivate the patient in the early days of therapy.
- Performance Check simplifies in-home device evaluation for easy remote troubleshooting.
- Connectivity: DreamStation offers more remote monitoring options than any leading CPAP brand.

DreamStation Auto BiPAP specifications

Part number	Description
INX700S15	DreamStation Auto BiPAP with SD card
INX700H15	DreamStation Auto BiPAP with humidifier and SD card

Fixed CPAP, fixed BiPAP, Auto BiPAP
EZ-Start
0 to 3
4-25 cm H ₂ O
0 to 45 min (5-min increments)
Patient adjustable
Heated humidification: fixed, adaptive, heated tube modes Heated tube (optional)
Color LCD/control wheel/push button
Ambient light sensor
Reusable pollen and disposable ultra-fine
Display: compliance 1 day, 7-day, 30-day averages SD card > 1 year On-board: up to 6 months
Automatic
2 years
15,7 (L) \times 19,3 (W) \times 8,4 cm (H) (without humidifier) 29,7 (L) \times 19,3 (W) \times 8,4 cm (H) (with humidifier)
1,33 kg (without humidifier), 1,98 kg (with humidifier) but including power supply
80W/12V/6,67A
100-240 VAC, 50/60 Hz, 2.0-1.0 A

Sleep therapy accessories	
Humidification	Part number
DreamStation heated humidifier	EUXH
Water chamber	1122520
Tubing	Part number
Heated tube (15 mm)	HT15
Standard tube (15 mm)	PR15
Filters	Part number
Pollen (foam) filter, reusable (1/pack)	1122446
Ultra-fine filter, disposable (1/pack)	1122447
Ultra-fine filter, disposable (2/pack)	1122518
Ultra-fine filter, disposable (6/pack)	1122519
Accessories	Part number
DreamStation carrying case	1121162
DreamStation link module/SpO ₂	1120293
DreamStation SpO ₂ cable	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, Adult, Flex	953
Flexiwrap, sens tape, Adult (25/pack)	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
DreamStation, 80W power supply	1118499
Connectivity	Part number
DreamStation cellular/SpO ₂ modem	100621C
DreamStation Wifi/SpO ₂ modem	100721W
· · · ·	





Philips Respironics 1010 Murry Ridge Lane Murrysville, PA 15668 USA

www.philips.com/respironics



DreamStation Auto CPAP features

- · Advanced event detection*
- Clinically validated auto algorithm (responds to obstructive events: AO, H, RE, FL, VS.)
- · Proactive search algorithm
- · CPAP-Check
- · Auto-Trial
- Opti-Start
- · NEW EZ-Start
- · NEW Smart Ramp (0-45 min)
- . Λ_ΕΙον
- · NEW Performance Check
- Encore patient management system: compliance and data reports, waveforms reports
- NEW Multiple remote monitoring solutions including cellular, Bluetooth® and Wifi modems
- · On-board memory
- · Oximetry module (optional)

Humidification

- $\boldsymbol{\cdot}$ Active humidification controlled by the room environment
- · Pre-heat humidifier capable
- · Heated tube (optional)

DreamStation Auto CPAP benefits

- **EZ-Start's** automatic, personalized adjustments to CPAP pressure is designed to help patients gradually acclimate to their prescribed level of therapy.
- Daily Progress Feedback displays the length of sleep and progress, designed to motivate the patient in the early days of therapy.
- Performance Check simplifies in-home device evaluation for easy remote troubleshooting.
- Connectivity: DreamStation offers more remote monitoring options than any leading CPAP brand.

DreamStation Auto CPAP specifications

Part number	Description
INX500S15	DreamStation Auto CPAP with SD card
INX500H15	DreamStation Auto CPAP with humidifier and SD card

Specifications	
Therapy modes	Auto, CPAP, CPAP-Check
Therapy features	Auto-Trial, EZ-Start, Opti-Start
Flex pressure relief	0 to 3
Pressure range	4-20 cm H ₂ O
Smart Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	Patient adjustable
Humidification	Heated humidification: fixed, adaptive, heated tube modes Heated tube (optional)
Device controls	Color LCD/control wheel/push button
LCD backlight control	Ambient light sensor
Filters	Reusable pollen and disposable ultra-fine
Data storage capacity (mininum)	Display: compliance 1 day, 7-day, 30-day averages SD card > 1 year On-board: up to 6 months
Altitude compensation	Automatic
Warranty	2 years
Dimensions	15,7 (L) \times 19,3 (W) \times 8,4 cm (H) (without humidifier) 29,7 (L) \times 19,3 (W) \times 8,4 cm (H) (with humidifier)
Weight	1,33 kg (without humidifier), 1,98kg (with humidifier) but including power supply
Power supply	80W/12V/6,67A
Electrical requirements	100-240 VAC, 50/60 Hz, 2.0-1.0 A


Sleep therapy accessories	
Humidification	Part number
DreamStation heated humidifier	EUXH
Water chamber	1122520
Tubing	Part number
Heated tube (15mm)	HT15
Standard tube (15mm)	PR15
Filters	Part number
Pollen (foam) filter, reusable (1/pack)	1122446
Ultra-fine filter, disposable (1/pack)	1122447
Ultra-fine filter, disposable (2/pack)	1122518
Ultra-fine filter, disposable (6/pack)	1122519
Accessories	Part number
DreamStation carrying case	1121162
DreamStation link module	1120293
DreamStation SpO ₂ module	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, Adult, Flex	953
Flexiwrap, sens tape, Adult (25/pack)	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
DreamStation, 80W power supply	1118499
	Part number
Connectivity	T di t Hullibei
Connectivity DreamStation cellular/SpO ₂ modem	100621C





Philips Respironics 1010 Murry Ridge Lane Murrysville, PA 15668 USA

www.philips.com/respironics



DreamStation helps support long-term patient use while creating **new efficiencies for your business**.

Connectivity options offer more monitoring choices than any leading CPAP brand with integrated Bluetooth® on every device, and offering optional Wi-Fi and cellular modems.*

Performance Check simplifies in-home device evaluation and can help reduce the number of normally functioning devices that are returned for service.

EZ-Start's automatic, personalized adjustments to CPAP pressure can help patients gradually acclimate to their prescribed level of therapy.

Daily Progress Feedback displays a simple trend of patients' nightly use to help track progress on the path to adherence.

DreamStation features

- · Advanced event detection
- · Mask type resistance control
- · Mask fit check
- · Oximetry capable
- Encore patient management system with optional Adherence Profiler
- $\cdot \ \mathsf{DreamMapper} \ \mathsf{patient} \ \mathsf{self-management} \ \mathsf{system} \ \mathsf{compatibility}$

Humidification

- · Fixed, adaptive, heated tube
- Pre-heat humidifier option
- · Modular design with one-piece, dishwasher-safe water chamber

Protection programs

- · Flex Promise program
- Bi-level Rescue program
- "One-Touch" warranty programs

BiPAP Pro features

- C-Flex, C-Flex+, Bi-Flex
- SmartRamp
- EZ-Start
- · CPAP-Check
- · Advanced event detection
- Waveform data
- · Standard integrated Bluetooth

*Internal assessment of 2015 competitive CPAP data comparing to ResMed Airsense10/Aircurve10 platform and Fisher & Paykel Icon series platform.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and the use of such marks by Philips Respironics is under license. Other trademarks and trade names are those of their respective owners.

DreamStation BiPAP Pro

	Comfort features		Encore data captur	e and reporting	
Modes	Features	Flex Comfort	Advanced Event Detection	Flow Waveforms	
Fixed CPAP Fixed BiPAP	EZ-Start	C-Flex, C-Flex+, Bi-Flex	Standard	Standard	

Specifications

DreamStation BiPAP Pro		
Pressure range	4 to 25 cm H2O	
Flex pressure relief	0 to 3	
Ramp time	0 to 45 min (5-minute increments)	
Starting ramp pressure	4 cm H2O to CPAP Pressure (CPAP); 4cm H2O to EPAP (BiPAP)	
Humidification	Heated humidification: Fixed, Adaptive, Heated Tube modes; Heated Tube (optional)	
Data storage capacity (minimum)	SD card: > 1 year On-board: <6 months	
Altitude compensation	Automatic	
Filters	Reusable pollen; disposable ultra-fine	
Device controls	LCD, control dial/push button	
Warranty	2 years (US)	
Dimensions	15.7x19.3x8.4cm (w/o humidifier) 29.7x19.3x8.4cm (w/ humidifier)	
Weight	1.33kg/2.94lbs (w/o humidifier) 1.98kg/4.37lbs (w/humidifier) Incl. power supply	
Electrical requirements	100 – 240 VAC, 50/60 Hz, 2.0-1.0 A	

Ordering information

DreamStation BiPAP Pro	Part number
DreamStation BiPAP Pro	INX600S15
DreamStation BiPAP Pro with humidifier	INX600H15
DreamStation BiPAP Pro with humidifier and heated tube	INX600T15
Power supply	Part number
DreamStation 80-watt power supply	1118499
Power cord	1005894

Ordering information

Humidification	Part number
DreamStation humidifier	DSXH
DreamStation water chamber	1122520
RP-DS Humidifier Dry Box Assembly	1120668
RP-DS Humidifier Flip Lid Seal	1120617
RP-DS Humidifier Dry Box Inlet Seal	1120613
Tubes	Part number
15MM Std Tube – DreamStation - RP	PR15
15MM Heated Tube	HT15
Filters	Part number
Pollen filter, reusable (1 per pack)	1122446
Ultra-fine filter, disposable (1 per pack)	1122447
Ultra-fine filter, disposable (2 per pack)	1122518
Ultra-fine filter, disposable (6 per pack)	1122519
Accessories	Part number
Link module	1120293
Nonin SpO2, assembly	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, adult, flex	953
Flexiwrap, sens tape, adult, 25/pk	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
Modems	Part number
Cellular modem, oximetry capable	100600C
Cellular modem, non-oximetry	100605C
Wi-Fi accessory, oximetry capable	100700W
Wi-Fi accessory, non-oximetry	100705W



^{© 2015} Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. (Royal Philips) or their respective owners.



DreamStation CPAP Pro features

- · Advanced event detection*
- · Therapy modes: CPAP, CPAP-Check
- · NEW EZ-Start
- · NEW Smart Ramp (0-45 min)
- · C-Flex
- · **NEW** Performance Check
- Encore patient management system: compliance and data reports, waveforms reports
- **NEW -** Multiple remote monitoring solutions including cellular, Bluetooth® and Wifi modems
- · On-board memory
- · Oximetry module (optional)

Humidification

- $\boldsymbol{\cdot}$ Active humidification controlled by the room environment
- · Pre-heat humidifier capable
- Heated tube optional

DreamStation CPAP Pro benefits

- **EZ-Start's** automatic, personalized adjustments to CPAP pressure is designed to help patients gradually acclimate to their prescribed level of therapy.
- Daily Progress Feedback displays the length of sleep and progress, designed to motivate the patient in the early days of therapy.
- Performance Check simplifies in-home device evaluation for easy remote troubleshooting.
- Connectivity: DreamStation offers more remote monitoring options than any leading CPAP brand.

DreamStation CPAP Pro specifications

Part number	Description
INX400S15	DreamStation CPAP Pro with SD card
INX400H15	DreamStation CPAP Pro with humidifier and SD card

Specifications	
Therapy modes	CPAP, CPAP-Check
Therapy features	EZ-Start
Flex pressure relief	0 to 3
Pressure range	4-20 cm H ₂ O
Smart Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	Patient adjustable
Humidification	Heated humidification: fixed, adaptive, heated tube modes Heated tube (optional)
Device controls	Color LCD/control wheel/push button
LCD backlight control	Ambient light sensor
Filters	Reusable pollen and disposable ultra-fine
Data storage capacity (mininum)	Display: compliance 1 day, 7-day, 30-day averages SD card > 1 year On-board: up to 6 months
Altitude compensation	Automatic
Warranty	2 years
Dimensions	15,7 (L) x 19,3 (W) x 8,4 cm (H) (without humidifier) 29,7 (L) x 19,3 (W) x 8,4 cm (H) (with humidifier)
Weight	1,33 kg (without humidifier), 1,98kg (with humidifier) but including power supply
Power supply	80W/12V/6,67A
Electrical	100-240 VAC, 50/60 Hz, 2.0-1.0 A

Sleep therapy accessories	
Humidification	Part number
DreamStation heated humidifier	EUXH
Water chamber	1122520
Tubing	Part number
Heated tube (15mm)	HT15
Standard tube (15mm)	PR15
Filters	Part number
Pollen (foam) filter, reusable (1/pack)	1122446
Ultra-fine filter, disposable (1/pack)	1122447
Ultra-fine filter, disposable (2/pack)	1122518
Ultra-fine filter, disposable (6/pack)	1122519
Accessories	Part number
DreamStation carrying case	1121162
DreamStation link module	1120293
DreamStation SpO ₂ module	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, Adult, Flex	953
Flexiwrap, sens tape, Adult (25/pack)	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
DreamStation, 80W power supply	1118499
Breamstation, cow power supply	
	Part number
Connectivity	
	Part number 100621C 100721W

© 2015 Koninklijke Philips N.V. All rights are reserved. Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time $\stackrel{,}{\underset{}{\stackrel{}}{}}$ without notice or obligation and will not be liable for any consequences resulting from the use of this publication.



Philips Respironics 1010 Murry Ridge Lane Murrysville, PA 15668 USA

www.philips.com/respironics