





# HbO-Smart

The new HbO-Smart device is the latest finger tip pulse oximeter from Drive DeVilbiss International intended for use in homes or hospitals for non-invasive measurement of oxygen saturation, pulse rate and perfusion index.

Including adjustable settings for user alerts and four optional screen displays this new stylish device has everything you would expect from a modern day spot checking oximeter which can be used for both children and adults.

## TECHNICAL SPECIFICATIONS

<b>Display Type</b>	OLED
<b>Dimensions</b>	65mm x 35mm x 30mm
<b>Oxygen Saturation</b>	
Range	70 - 99%
Accuracy	80 - 99% / ±2%
Resolution	1%
<b>Pulse Rate</b>	
Range	30 bpm - 240 bpm
Accuracy	± 2 %
Resolution	1 bpm
<b>Perfusion Index</b>	
	0.3 - 20%
<b>Battery</b>	2 x AAA
<b>Battery Life</b>	25 hours
<b>Power Consumption</b>	< 30mA



SpO<sub>2</sub>, Pulse Rate and Pulse Rate Bar, Plethysmograph Waveform, Perfusion Index and Battery Status Indicator

The HbO-Smart is ready to use as it includes a lanyard, a 5 language user manual and 2 x AAA batteries.



### Drive DeVilbiss International

DeVilbiss Healthcare GmbH  
Kamenzer Straße 3  
68309 Mannheim, Germany

UK +44 (0) 121 / 521 31 40  
DE +49 (0) 621 / 1 78 98 0  
kontakt@devilbisshc.com

[www.drivedevilbiss-int.com](http://www.drivedevilbiss-int.com)

