



Fit more patients comfortably

Respironics PerformaTrak Full Face Mask for single-patient use (EE, SE)

The Respironics PerformaTrak Full Face Mask minimizes facial pressure while maintaining a good seal, even at the nasal bridge. This helps resolve problems associated with the application of noninvasive ventilation. By easily adjusting the mask force at the bridge of the patient's nose, clinicians can provide a patient with the comfortable fit needed to enhance ventilation therapy.

Key advantages

- Dual-density foam forehead pads provide patient comfort and mask stability.
- With only three sizes, clinicians can quickly fit a wide range of patients.
- Swivel clips simplify re-application by eliminating the need to refit after removal.

PHILIPS

Respironics PerformaTrak

The Respironics PerformaTrak Full Face Mask has additional advantages to improve patient comfort and noninvasive ventilation (NIV) performance.

- Adjustable forehead support improves overall mask fit and reduces skin breakdown.
- With a choice of the standard elbow (SE) for dual-limb systems and the entrainment elbow (EE) for single-limb systems, clinicians can use the PerformaTrak Mask with any noninvasive ventilator.
- Digital Auto-Trak compatibility enhances NIV performance.
- Integrated sizing gauge provides a convenient way to reduce cross-contamination and waste.
- Low cost replacement headgear is available in two sizes
- Velcro® grab tabs allow quick, easy removal and adjustment of headgear.
- Optional sealing pad minimizes air leaks and localized skin breakdown from nasogastric (NG) tubes.

Description	Quantity	Part Number
PerformaTrak Full Face Mask, small	1	1012623
PerformaTrak Full Face Mask, medium	1	1012624
PerformaTrak Full Face Mask, large	1	1012572
PerformaTrak Full Face Mask, small	10	1018592
PerformaTrak Full Face Mask, medium	10	1018593
PerformaTrak Full Face Mask, large	10	1018594
PerformaTrak head strap with clips	1	1015788
Nasogastric tube sealing pad	10	452040



Entrainment elbow (EE)



Standard elbow (SE)



© 2009 Koninklijke Philips Electronics N.V. All rights are reserved. PerformaTrak and Respironics are trademarks of Respironics, Inc. and its affiliates.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com
 fax: +31 40 27 64 887

Hoech JP MCI 4102185 PN 1060155
 4522 962 44101 * FEB 2009

Philips Healthcare
 Global Information Center
 P.O. Box 1286
 5602 BG Eindhoven
 The Netherlands