



Experience the Accuracy

MULTI PARAMETER MONITOR

EXCELLO / EXCELLO M / EXCELLO D



Happier Living Everyday



The device provides arrhythmia detection, ST & HRV analysis with 3-lead/5-lead ECG monitoring and comprehensive trend data review



High resolution Display

High resolution 10.4" colour TFT display



Standard Configuration

Standard configuration of ECG, RESP, Nellcor SpO₂ / Masimo SpO₂ / Digital SpO₂, NiBP, TEMP



Drug Dose Calculation

The added functions of the device are drug dose calculation / HRV analysis / OXY-CRG



Arrhythmia and ST Analysis

The device provides for arrhythmia detection and ST analysis



ECG Monitoring

5 Lead ECG / 3 Lead ECG (Optional)



Waveform Display

Maximum 10 waveforms display



168 Hour Trend

Provides for maximum of 168 hours of trend review



Pacemaker Detection

Detects pacemaker pulses



Rechargeable Battery

>120 minutes operating time in fully charged condition and >240 minutes with additional battery (optional)



Networking

Wired or Wireless networking with Central monitoring system software (optional)



Optional Parameters & Configurations

2-IBP, ETCO₂ (mainstream/side stream/ micro stream), AGM, ICG and thermal recorder

*CE marking is applicable for EXCELLO

CERTIFIED ISO 13485:2003, ISO 9001:2008 COMPANY

BPL Medical Technologies Private Limited
11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076
CIN: U33110KA2012PTC067282
Ph: +91 80 26484388/ 2648 4350
Toll Free: 1800-4252355
Website: www.bplmedicaltechnologies.com
For Enquiries: sales.medical@bpl.in



A BPL Medical Technologies Magazine



Follow Us on

