



PORTABLE OXYGEN CONCENTRATOR


Product Code : MKPOC01



Battery Duration

Setting flow	1 : 5 h 24 min
	2 : 4 h 48 min
	3 : 3 h 18 min
Battery duration	4 : 2 h 30 min
	5 : 2 h 06 min
	6 : 1 h 42 min

Powerful Compressor



Specifications :

Oxygen Concentration	90%-3%/+6% at all settings
Oxygen Flow Max	1.2 Liter/min
Warm-Up Time	2 minutes
Weight	4.37 lbs. 1.98kg (includes battery)
Size	Length: 8.70" (221mm) Width: 6.30" (160mm) Height: 3.35" (85mm) (includes battery)
Oxygen Flow	Pulse dose delivery system Five flow settings 1, 2, 3, 4, 5, 6
Power	Battery: Voltage: 14.4VDC Rated capacity: 6.7Ah AC Input: 100to240VAC 50 to 60 Hz DC Input: 11-16VDC DC Output: 19V 6.3A Duration : Up to 5.4 hours
Battery	Recharge approximately 4 hours with AC or DC power
Noise Level	49 Decibels* @setting 2
Display Screen	2.8 inch
Warranty	3years for main unit ,1 year for Sieve Bed and Battery
Operation	Simple control functions and easy-to-read LCD display
Use	Designed for 24/7 use



High Efficiency



Excellent Heat Dissipation



Low Noise



Long Life

