

Oxy9Wave


The Most Accurate Oximeter

High Performance

- * 3.2" Color TFT LCD for Clear View
- * Large memory : Up to 30 Days in 10 Secs
- * Nellcor compatibility

User Convenience

- * Vertical & horizontal display of measurement value
- * Intuitive operation



Vertical display



Horizontal display



SPECIFICATIONS

<ul style="list-style-type: none"> ■ Performance 	<p>SpO₂ (Blood oxygen saturation)</p> <p>Range : 40 to 100% Accuracy : 70~100%, ±2% 0~69% Unspecified</p> <p>BPM (Heartbeat)</p> <p>Range : 30 to 250 bpm Accuracy : ±2bpm</p> <p>Trend</p> <p>Memory : 30 days at 10 second interval</p>
<ul style="list-style-type: none"> ■ Display 	<p>3.2" TFT LCD (Resolution : 240 × 400 pixel)</p>
<p>Display parameter</p>	<p>SpO₂, BPM, SpO₂ limit, BPM limit, Level Bar, Heartbeat wave, Time and Date, Alarm volume, Battery status and brightness</p>
<p>LED Display</p>	<p>Sensor separation, Alarm silence, Charging status and Power connection</p>
<ul style="list-style-type: none"> ■ Power Supply 	<p>24VA (max)</p>
<p>Adaptor</p>	<p>Input : 100-240V~, 50-60Hz, 0.4A Output : +5V, 2.0A</p>
<p>Battery</p>	<p>Type : Lithium-ion battery 4400mAh, 3.7V Operating time : Up to 20 hours Charging time : 6.5 hours</p>
<ul style="list-style-type: none"> ■ General 	<p>Operating temperature</p> <p>5°C to 40°C</p> <p>Non-operating temperature</p> <p>-20°C to 70°C</p> <p>Humidity</p> <p>5% to 95%</p> <p>Operating altitude</p> <p>-1,000 ft to 12,000 ft</p> <p>Weight</p> <p>300g</p> <p>Size</p> <p>45(H) × 210(W) × 75(D)mm</p>
<ul style="list-style-type: none"> ■ Accessory 	<p>Adult reusable SpO₂ finger sensor 1ea AC adaptor 1ea Pole mount bracket 1ea</p>
<p>Standard</p>	<p>Neonate reusable SpO₂ finger sensor (Y-type) Extension cable</p>
<p>Option</p>	



Main office : Bionet Co.,Ltd.

Tel : +82-70-7585-6418 (Int'l)

e-mail : sales@ebionet.com (Int'l)

North America : Bionet America, Inc.

Tel : +1-714-734-1760 (U.S.A.)

Europe : MGB Endoskopische Geräte GmbH Berlin

Tel: +49-30-6392-7000 (Europe)

www.ebionet.com

02-6300-6477 (대한민국)

yunms@ebionet.com (대한민국)

www.bionetus.com

e-mail:sales@bionetus.com (U.S.A.)

www.mgb-berlin.de

e-mail:info@mgb-berlin.de (Europe)