



ResMed



13.6 cm, real size

AirMini™ Solutions

Rewrite your bedtime story

ResMed.com/AirMini

AirMini: effective¹⁻⁴ compact therapy without compromise

Made by ResMed, you know that sleep therapy on AirMini Solutions will be just as good as the therapy delivered by our other market-leading devices - in a compact, friendly, easy-to-use and more appealing package.

AirMini Solutions provide effective¹⁻⁴ therapy at their most comfortable and mobile. You don't have to compromise between therapy performance, better health and freedom anymore.




HumidX™ waterless humidifier

Compact and discreet, HumidX* or HumidX Plus* is the ideal choice for patients looking to experience therapy comfort on the go and not only.




ActiveAir™

This crucial technology regulates airflow to keep pressure delivery stable and comfortable at all levels.



* HumidX and HumidX Plus are only compatible with AirFit™ N20 and AirFit P10 for AirMini.



AirFit™ and AirTouch™ series of masks

AirMini is compatible with N20 and F20 series of masks, and the AirFit P10 for AirMini.

Start & Stop


With only one Start button, it's so easy to use your device every night.

Bluetooth connection

Just press the Bluetooth button, and synchronise your AirMini with your smart device.

AirMini App

A smart device app that allows you to manage your comfort settings and view your therapy progress. The AirMini App helps ensure you are swiftly set up for therapy.



Power to go, freedom to move

Because AirMini Solutions can deliver effective¹⁻⁴ therapy without the bulk, you can take it from the bedside table to the hotel room - and even to the great outdoors. It's simply small, convenient and smart.



Freedom of choice

AirMini Solutions are available with a range of accessories and customisation options so that therapy can be personalised to suit various individuals' different needs and preferences.

These options include various algorithms, therapy modes, features, as well as mask shapes and styles and accessories.



Small

AirMini is the smallest⁵ CPAP device ever manufactured.

Smart

The AirMini App allows you to customise your therapy comfort to your individual needs and view your progress with ease.

Proven¹⁻⁴

AirMini features the same established algorithms as ResMed's AirSense range of sleep therapy devices.

Trusted¹⁻⁴ therapy made compact

Understanding your everyday needs

When it comes to your health, you shouldn't have to compromise on anything: comfort, ease-of-use, effective¹⁻⁴ therapy, humidification, performance, size and mobility. These are all the standard features you have the right to expect from your therapy, and that AirMini delivers.



Technology that serves you

We have developed a unique venting technology to deliver market-leading therapy, portable waterless humidification and uncompromised performance⁶ to you, all in a small/light and lifestyle-friendly package.

Engage with your well being with the AirMini App

Track and manage your own sleep quality

With AirMini Solutions, you have access to ResMed's AirMini patient engagement app. Whether you're a new or trained user, the AirMini App allows you to track your sleep data and manage your therapy from the convenience of your smart device⁶. Very friendly and quick, with easy connection to your machine.

Get support at every stage of therapy

The AirMini App also provides coaching and timely personalised support, and gives you access to self-help tools to help you manage your therapy, care for your device and ease into treatment.



Download app
Install your AirMini App on your smart device⁶



Sync your AirMini
Press the Bluetooth button and activate the bluetooth option on your smart device



Adjust features
Adjust comfort features from your smart device (Ramp or EPR)



Track your sleep
Check out your sleep score every morning



Tips & coaching
Learn ways to improve your comfort and sleep better



Download the AirMini App now!
...and stay connected

Keeping you and your healthcare team connected, your detailed therapy data can be sent to your clinician/healthcare provider to review as part of your ongoing therapy management and care.





Travel accessories to make your journey easier

With a variety of optional accessories, the AirMini system is designed to suit your lifestyle.

A hard case offers protection – while a travel bag allows you to transport the entire AirMini system (tubing, mask and device) with ease. The DC/DC Converter 65W allows you to operate your machine from a vehicle cigarette lighter socket.

The AirMini mount system even allows your AirMini to be secured in a variety of bedside configurations for greater convenience and discretion.

Visit mysleep.resmed.com to discover all these accessories.



visit ResMed.com/AirMini

or buy online
mySleep.resmed.com

- 1 Netzel et al. APAP device technology and correlation with patient compliance. *Somnologie - Schlafforschung und Schlafmedizin* 2014; DOI: 10.1007/s11818-014-0662-0.
- 2 Isetta et al. Comparative assessment of several automatic CPAP devices responses: a bench test study. *ERJ Open Res* 2015; 1:00031-2015.
- 3 Zhu et al. All APAPs are not equivalent for the treatment of sleep disordered breathing: a bench evaluation of eleven commercially available devices. *J Clin Sleep Med* 2015 11(7):725-34.
- 4 Isetta et al. Novel Approach to Simulate Sleep Apnea Patients for Evaluation Positive Pressure Therapy Devices. *PLoS ONE* 11(3): e0151530.doi:10.1371/journal.pone.0151530.
- 5 As of 04, April 2017. 136 mm (W) x 84 mm (H) x 52 mm (D), 250 cm³.
- 6 Equivalent to mid-range AirSenseTM HumidAir setting. Internal benchtesting – data on file