








# Biocard® HCV RAPID TEST KIT


Biocard® HCV Rapid is qualitative in-vitro diagnostic test for the detection of antibodies to Hepatitis C Virus

-  Approved by National Institute Biologicals (NIB)
-  Sensitivity : 100% Specificity : 100%
-  Presence of control line to validate the assay
-  Sample types: Serum, Plasma (10µL)/ Whole Blood (20µL)
-  Pre coated with HCV Core NS3, NS4 and NS5
-  Results with in 20 minutes
-  Storage: 4°C to 30°C




# Biocard® HCV RAPID TEST KIT

## Test Procedure



## Results Interpretation

For Biocard® HCV



\*refer kit insert for exact procedure

## Performance

Status	Positive	Negative
Positive	110	0
Negative	0	650
Total	110	650

Diagnostic Sensitivity: 100 %

Diagnostic Specificity: 100 %

## Ordering Information

Product Code	Product Name	Pack Size	Storage	Shelf Life
30201416	<b>Biocard®</b> HCV RAPID TEST KIT	30 Test	4-30°C	24 Months

\*Please note: Product specifications are subject to change without prior notice owing to product modifications, improvements / up-gradation.



corporate@trivitron.com | www.trivitron.com

Customer Care Helpline: +91 98400 80008



**TRIVITRON**  
HEALTH CARE  
speaking your language