










Biocard[®] SYPHILIS RAPID TEST KIT


Biocard[®] Syphilis Rapid test kit is an immunochromatographic immunoassay for the qualitative detection of antibodies (IgG, IgM and IgA) to *Treponema pallidum*

-  Approved by National Institute Biologicals (NIB)
-  Sensitivity : 100% Specificity : 100%
-  Presence of control line to validate the assay
-  Sample types: Serum, Plasma (10µL)/ Whole Blood (20µL)
-  Pre-coated with recombinant Treponemal antigens (Tp15, Tp17, and Tp47)
-  Results with in 20 minutes
-  Storage: 4°C to 30°C




Biocard® SYPHILIS RAPID TEST KIT

Test Procedure



Results Interpretation

For Biocard® Syphilis



Performance

Status	Positive	Negative
Positive	75	0
Negative	0	580
Total	75	580

Diagnostic Sensitivity: 100 %

Diagnostic Specificity: 100 %

Ordering Information

Product Code	Product Name	Pack Size	Storage	Shelf Life
30201417	Biocard® SYPHILIS RAPID TEST KIT	30 Test	4-30°C	24 Months

*Please note: Product specifications are subject to change without prior notice owing to product modifications, improvements / up-gradation.



corporate@triviron.com | www.triviron.com

Customer Care Helpline: +91 98400 80008



TRIVITRON
HEALTH CARE
speaking your language