



Coag4D

4-CHANNEL COAGULOMETER

- Compact design & versatile 4-channel Semi-automated coagulometer
- Complete test palette for clotting, immunoturbidimetric and chromogenic tests
- Supports D-Dimer and AT-III tests
- Easy-to-use operation, Innovative user interface with built-in QC feature



SPECIFICATIONS

Key Features:

- 4 independent measuring channels (2 parallels or 4 single tests)
- Screening tests, Factors, Inhibitors & User defined tests
- Quantitative D-Dimer measurement
- AT-III activity measurement
- Maintenance free
- Liquid crystal color touch screen
- Innovative user interface
- Visual sample traffic guidance for easy and touchless operation
- Synchronized on screen and LED real-time channel status follow-up
- QC feature
- On-board storage of up to 20 cuvettes for preincubation

- Green for preparation
- Orange for incubation
- Red for measurement



System Elements:

- 2 thermostated reagent positions both with magnetic stirring
- Automatic Thermal Printer, 2 USB connections for software upgrade and data transfer
- Optional barcode reader
- 405 nm, 570 nm and 640 nm wavelength detection

*Please note: Product specifications are subject to change without prior notice owing to product modifications, improvements / up-gradation.