



## Coaxial Ophthalmoscope Rechargeable Set

- Facilitates ophthalmoscopy by enabling easier entry into the eye, a larger field of view and reduced glare compared to standard ophthalmoscopes
- LED (SureColour) light for true tissue colour and consistent long-lasting performance
- Produce a shadow-free spotlight, easier entry into undilated pupils and a larger field of view
- Equipped with 28 focusing lenses with a range of -25 to +40 diopters
- It has cross-linear polarising, red-free and cobalt blue filters
- It has micro-spot, small spot, large spot, slit and fixation aperture
- Features a smooth, non-abrasive knurled finished rechargeable Nickel-Cadmium 60 minute power handle that can also be converted for use with standard C-cell batteries

**Ordering Information:** 11772-VCL-3.5V Coaxial Ophthalmoscope Rechargeable Set