

MHL 1008-C aesthetically designed ultralight aluminum alloy wheelchair with strong, flexible, fire-retardant and antibacterial upholstery providing comfort of 17.7 inches' width.



MHL 1008-C

LIGHTWEIGHT MANUAL ALUMINIUM WHEELCHAIR

MHL 1008-C

LIGHTWEIGHT MANUAL
ALUMINIUM WHEELCHAIR

DIMENSIONS

Overall Length	905 mm
Overall Height	870 mm
Overall Width Folded	260 mm
Overall Width Open	600 mm
Seat Height From Footrest	440 mm
Seat Height From Ground	480 mm
Armrest Height From Seat	228 mm
Ground Clearance	145 mm
Ground Clearance from Footrest	95 mm
Front Wheel	6"
Rear Wheel	16"x1" 3/8
Seat Width	450 mm

GENERAL SPECIFICATIONS

Propulsion System	Manual
Frame	Aluminium frame with MS cross tubes
Upholstery	Canvas, Antimicrobial & Fire-retardant
Weight Capacity	Up to 100 Kg
Net Weight	9.8 Kg
Front Wheel	Mag Design PP Castor
Rear Wheel	MAG Wheel with SOLID TYRE
Foldable / Fixed	Foldable
Foot Rest	Polymer And Flipable
Back Rest	Flip/Foldable
Armrest	PU-Moulded, Non-detachable
Head Rest	No
Anti Tip Provision	No
Propelling Ring	Yes with plastic
Manual Brakes	Anti Skid Braking System
Attendant Brakes	Yes
Push Handles	Yes with PVC Grip
Provision Of Commode	No
Calf Strap	Yes



Cushioned arm rest gives ease of handling and comfort. Smooth powder coated propelling ring for strong grip.



Foldable backrest, calf support and detachable foot rests makes it perfect for senior citizens and people with compromised mobility, and folds within 11.2 inches of width.



Rear wheels of 16 inches for ease of maneuverability and solid PP front wheels of 6 inches with double cross design enables 100 kg of weight bearing strength. Small storage space in the form of back – pocket for storage of documents and medicines.



HEROeco
MED

#57, Udyog Vihar Phase IV, Sector 18, Gurugram 122015,
Haryana, India | Tel: +91 124 683 0000