



HEROeco
MED

MHL 4001 is an ideal solution of transport for specially-abled. This is a double-hand drive tricycle that easy to handle. The frame and footrest is made with sturdy MS Steel that makes it very durable and long-lasting.



MHL 4001
TRICYCLE

MHL 4001

TRICYCLE

DIMENSIONS

Overall Length	1960 \pm 50 mm
Overall width	890 \pm 25 mm
Overall height	990 \pm 25 mm
Width of footrest	320 \pm 10 mm
Length of footrest	600 \pm 10 mm
Clearance of footrest from ground	145 \pm 25 mm
Armrest height from seat	225 \pm 10mm
Seat length	600 \pm 10 mm
Seat width	430 \pm 10 mm
Back height from seat	320 \pm 10 mm
Length of steering handle	420 \pm 10mm
Height of rear wheel supporting frame	385 \pm 10mm
Leverage of steering handle	5:01



It comes with a steering handle to maneuver the tricycle giving more control over different road conditions and makes it easy to ride. The footrest is 12.6 inches wide giving ample space for the user.

Heavy-duty wheel rims, mudguard, driving chains are used to make this tricycle strong and durable. It comprise of two hand powered steerable powered rear wheels and one coasting front wheel. The crank mechanism makes it easy to power the tricycle using both hands when the hand levers are rotated.

The seat is well-cushioned and fitted perfectly on the frame to give a comfortable experience to the user. The tricycle is anticorrosive and comes a standard black color



#57, Udyog Vihar Phase IV, Sector 18, Gurugram 122015,
Haryana, India | Tel: +91 124 683 0000