



REHABILITATION AIDS

ABOUT US

We are a healthcare equipment manufacturer delivering high quality medical equipment for individual and institutional requirements. Our product design and equipment manufacturing meets Global Standards which caters to the growing Indian market in addition to worldwide exports.



MHL 1002 R – EVA / PU MANUAL CHROME PLATED
MAG WHEELCHAIR

MHL 1002 R – EVA/PU

MANUAL CHROME PLATED MAG WHEELCHAIR

DIMENSIONS

	EVA / PU
Overall Length	1045 mm
Overall Height	870 mm
Overall Width Folded	285 mm / 290 mm
Overall Width Open	685 mm
Seat Height from Footrest	460 mm
Seat Height from Ground	485 mm
Armrest Height from Seat	275 mm
Ground Clearance	70 mm
Ground Clearance from Footrest	75 mm
Front Wheel	200 mm
Rear Wheel	605 mm
Rear Wheel Tyre Width	36 mm
Seat Width	500 mm

GENERAL SPECIFICATIONS

Propulsion System	Manual
Frame	Chrome Plated Strong Frame
Upholstorey	Pvc , Antimicrobial
Weight Capacity	Upto 100 Kg
Net Weight	14.28 Kg / 14.40 Kg
Front Wheel	Mag Design PP Castor
Rear Wheel	MAG Wheel - Black EVA Tyre / MAG Wheel - Grey PU Tyre
Foldable / Fixed	Foldable
Foot Rest	Metallic and Detachable
Back Rest	Non Reclining / Foldable
Armrest	Cushioned, Non-detachable
Head Rest	NO
Anti Tip Provision	NO
Propelling Ring	Yes , Polymer
Manual Brakes	Anti Skid Braking System
Attendant Brakes	No
Push Handles	Yes With PVC Grip

HL 1002 R – EVA / PU is ergonomically designed chrome-plated wheelchair with strong, flexible and antibacterial upholstery providing comfort of 19 inches' width.



Foldable feature enhances the portability and can be stored easily in house or car.

MAG Rear wheels are available in two options of BLACK – EVA tyres and GREY – PU tyres of 24 inches for ease of maneuverability and solid PP front wheels of 8 inches enables 100 kg of weight bearing strength.



Smooth polymer-propelling ring for sound grip and reduced effort, optimized for 11'o clock to 2'o clock pushing ark.

Calf support and detachable metallic foot rests makes it optimal solution for patients with compromised mobility, neuromuscular issues, musculoskeletal problems and having fractured bones.



57, Udyog Vihar Phase IV, Sector 18, Gurugram 122015,
Haryana, India | **Tel.:** +91 124 683 0000