



Features :

- ▶ Portable delicate design , big Color Display, very easy to operate
- ▶ Accurate pulse pace identification function and excellent auto-interpretation function
- ▶ High accuracy digital filter, auto-baseline adjustment
- ▶ 12 lead synchronously acquisition and record
- ▶ Working modes : auto/Man/R-R Rhythm/Storage
- ▶ 800*480 7 inch LCD display detailed ECG information
- ▶ 250 cases storage and replay
- ▶ Rechargeable 2000mAh Li-ion battery supports continuous works for about 3-4 hrs
- ▶ Recording Mode : 3Ch, 3Ch+, 3Ch+++,6Ch
- ▶ Real-time freeze function
- ▶ 112mm*20m high resolution roll paper to print ECG
- ▶ RS232 and USB interface support USB disk storage ,

Optional :

PC ECG paper

Amigo Medical Systems

38/2953-6, 1st Floor, Joseph tower, CBI Road Kathrikadavu, Kaloor, Cochin- 682017

M: +91 8606141100 | +91 8606751100 | +91 9037995552

E: sales@amigomedicals.in | www.amigomedicals.in

