



## PulmoAide

### The perfect in-home asthma treatment

#### TECHNICAL SPECIFICATIONS

Dimensions (H x W x D)	28 cm x 27 cm x 17 cm
Weight	3.2 kg
Electrical Requirements	230 V / 50 Hz / 100 Watts
Maximum Pressure	2.1 ba or greater
Compressor Free Air Flow	9 l/min. or greater
Pressure Rate Nebuliser	0.8 - 1.1 bar
Nebuliser Litre Flow	5 - 7 l/min.
Nebuliser Capacity	6 ml (cc)
Nebuliser Rate	0.15 ml/min.
Particle Size / MMAD	75 % < 5 µm / 3.6 - 4.0 µm

The PulmoAide Nebuliser high pressure compressor from Drive DeVilbiss has been providing effective aerosol treatments for many years. Modern styling, durability and ease of use makes this unit perfect for in-home inhalation therapy.

The PulmoAide provides an unbeatable combination of quality and value designed to keep pace with today's competitive world of managed health care. Combined with the VixOne nebulizer, the PulmoAide's design and reliability make the performance breath enhanced Nebuliser the ideal choice for the hospital or home. You can trust.

The PulmoAide comes ready to use with all the attachments – filter, tubing, medicine cup and mouthpiece.

#### Product Features:

- Portable, lightweight and robust design
- Built-in storage compartment to hold nebuliser and tubing until needed again
- Integrated handle for easy carrying
- Easy set-up within seconds the PulmoAide is ready with an MMAD of 3.6 - 4.0 µm.
- High Aerosol output to reduce treatment time
- Spill proof design
- VixOne small volume jet nebulizer included
- Nebuliser can be cleaned in the dishwasher
- Mouthpiece can be replaced with a standard mask
- Less medication waste due to effective treatment and a residual volume of < 0.5 ml
- Gives effective treatment and results, even at an angle of 45°



Filter, tubing, nebulizer and mouthpiece included

#### Drive DeVilbiss International

DeVilbiss Healthcare GmbH  
Kamenzer Straße 3  
68309 Mannheim, Germany

UK +44 (0) 121 / 521 31 40  
DE +49 (0) 621 / 178 98 0  
kontakt@devilbisshc.com



www.drivedevilbiss-int.com

