

Ultrasound

1 & 3 MHz Ultrasound

JUS 2



Therapeutic Ultrasound for Every Clinic

Johari's ultrasound devices provide effective Continuous and Pulsed treatments in pain relief and muscle rehabilitation. It delivers 1 & 3 MHz through a unique ultrasound applicator. The effect of ultrasound increase in local blood flow can be used to help reduce local swelling, chronic inflammation and promote bone fracture healing.

Features

- Digital display
- 3.8 cm dia Ultrasound applicator
- Selectable Duty Cycle and Frequency
- Microprocessor controlled digital unit
- Inbuilt Intensity Control in applicator

Accessories of JUS-2



Ultrasound Applicator

Technical Specification

Outputs	- One	Max. Intensity	- 3 watt/cm ² (in step of 0.1 watt)
Waveform	- Ultrasound (Sinewave)	Treatment Time	- 1 - 60 Minutes (in step of 1 min)
Modes :	- 2 (Pulsed and Continuous)	Power	- 110 - 220 V AC
Frequency	- 1 MHz and 3MHz		
% Duty Cycle	- 10%, 20% and 50%		