

Romo Soft Kit™

[Home](#) » [Domestic](#) » [Products](#) » [Critical Anesthesia Care](#) » [Romo Soft Kit™](#)

- **Product Name:** Romo Soft Kit™
- **Generic:** PERCUTANEOUS TRACHEOSTOMY KIT
- **Product Code:** GS-2103
- **Product Category:** Critical Anesthesia Care



Romo Soft™ Kit (GS-2103)

- ▶ Specially designed kit permits a safe, fast, and cost-effective percutaneous tracheostomy procedure with minimal complications.
- ▶ Available with or without Dilating Forceps

COMPLETE KIT CONSISTS OF :

- ▶ Romo soft percutaneous tracheostomy tube
- ▶ Dilating forceps
- ▶ J-tip guide wire
- ▶ Syrin(ge-10) ml
- ▶ Dilating cannula
- ▶ Dilator
- ▶ Scalpel
- ▶ Cotton tape

SPECIAL FEATURES

- ▶ Bedside procedure
- ▶ Over the guidewire technique
- ▶ Forceps allow rapid dilation
- ▶ Accurate Placement

Tube I.D

7.0 mm

7.5 mm

8.0 mm

8.5 mm

Tube O.D

10.50 mm

11.30 mm

11.90 mm

12.60 mm

Box Quantity

1

Master Box Quantity

50