

# NewLife® Intensity® 10

## High Flow Stationary Concentrator



### Built-to-last Oxygen Delivery Innovation

Since being introduced in 2009, the New Life Intensity 10 has been a trusted high flow stationary concentrator addressing patients' oxygen delivery needs in over 100 countries including specialty applications in home and medical care facilities.

#### Time-proven design

- **Industry-leading reliability**—greater than 99% reliability rate out of the box and less than 2%\* warranty repair rate
- **Comprehensive warranty**—2-year warranty on the device
- **Durability**—engineered with an oversized compressor and sieve beds to minimize service

#### Versatile oxygen delivery

- **Oxygen delivery**—up to 95.5% concentration at up to 10 LPM, meeting a wide range of patient prescriptions
- **High outlet pressure**—the first oxygen concentrator offering 20 psi to support specialty applications and up to 200 feet of oxygen tubing
- **Energy efficient**—reduced power consumption, minimizing the burden on the patient's energy bill

#### Ease of service

- **Minimal maintenance**—proven quality design and performance to minimize required service
- **Modular design**—fewer components reduce time and effort to service the device
- **Enhanced diagnostics**—onboard alarm alerts enable quick troubleshooting

- ✓ US FDA Certified
- ✓ European Certificate –CE
- ✓ Indian CDSCO Certified

### NewLife Intensity 10 Specifications

**FLOW RATE**  
2.0 – 10.0 LPM

**POWER**  
600 W

**DIMENSIONS**  
27.5 x 16.5 x 14.5 inches  
(69.9 x 41.9 x 36.8 cm)

**MAXIMUM TUBING**  
200 ft (61 m)

**WEIGHT**  
58 lbs (26.3 kg)

**WARRANTY**  
2 years

**SOUND LEVEL**  
58 dB(A)

**MAINTENANCE SCHEDULE**  
Felt Filter – 1 Year Replacement  
Intake Filter – Clean Weekly

**OUTLET PRESSURE**  
20 psig (138 kPa)



Caire Inc.  
2200, Airport Industrial Dr. Ste.  
500 Ball Ground, GA 30107 USA  
<https://www.caireinc.com/product>