

VisionAire™ 5

Stationary Oxygen Concentrator



The foundation of your patients' oxygen therapy

Compact and lightweight oxygen delivery

The VisionAire 5 has been completely redesigned to be state-of-the-art, meeting all current compliance standards with quality engineering and reliability.

Time-proven design

- **Industry-leading reliability**—greater than 99%* reliability rate out of the box
- **Comprehensive warranty**—3-year warranty on the device
- **Proven durability**—engineered with robust valves and sieve beds to minimize failure rates

Promotes patient satisfaction

- **Effective oxygen delivery**—up to 95.5% concentration from 1.0 to 5.0 LPM, meeting a wide variety of patient prescriptions
- **Quiet**—you will hardly know this lightweight and compact device is there
- **Energy efficient**—reduced power consumption to minimize the burden on the patient's energy bill

Ease of service

- **Maintenance free**—no required preventative maintenance reducing your total cost of ownership
- **Dual hour meter**—clear LCD provides time of use in hour and minutes
- **Enhanced diagnostics**—on board alarm alerts facilitate the ability to quickly troubleshoot

VisionAire 5 Specifications

FLOW RATE
1.0 – 5.0 LPM

POWER
290 W (all flow rates)

DIMENSIONS
52.8 x 35.8 x 29.2 cm

MAXIMUM TUBING
15.2 m (50 ft)

WEIGHT
13.6 kg

WARRANTY
3 years

SOUND LEVEL
39 dB(A)

- ✓ US FDA Certified
- ✓ European Certificate –CE
- ✓ Indian CDSCO Certified

HUMIDIFIER BOTTLE

Item Number HU003-1



HUMIDIFIER ADAPTOR

(9 in. / 23 cm.)

Item Number 20843882



Caire Inc.
2200, Airport Industrial Dr. Ste.
500 Ball Ground, GA 30107 USA
<https://www.caireinc.com/product>