

# FreeStyle® Comfort®

The perfect fit



## Key features and benefits



### Optimal Comfort

Let your patients experience the comfort of the body-skimming, curved silhouette design. The FreeStyle Comfort portable oxygen concentrator is lightweight and compact.



### Smart O2 delivery

Unique clinical features including UltraSense®, autoDOSE®, and myCAIRE™ help to ensure your patients will get the oxygen they need, when they need it.



### Proper O2 Saturation

With 5 settings and up to 1050 mL of oxygen per minute, the FreeStyle Comfort serves a wide range of patient needs. The unit delivers up to 95% pure oxygen to further ensure proper saturation.



### Proven Quality

Robust design features, including the sieve bed technology and critical components, allow for minimum maintenance and long life. CAIRE offers one of the best warranties available, with a 2-year warranty on the device and sieve beds.

- ✓ US FDA Certified
- ✓ European Certificate –CE
- ✓ Indian CDSCO Certified
- ✓ FAA Approved –Airline Travel

## Specifications



### FLOW RATE

Pulse Dose 1-5  
O2 Output Settings:  
Setting 1–210 mL/min  
Setting 2–420 mL/min  
Setting 3–630 mL/min  
Setting 4–840 mL/min  
Setting 5–1050 mL/min

### DIMENSIONS

10.0 x 7.3 x 3.1 in  
(25.4 x 18.5 x 7.9 cm)  
3-cell battery pack  
11.0 x 7.3 x 3.1 in  
(27.9 x 18.5 x 7.9 cm)  
16-cell battery pack

### WEIGHT

5 lbs (2.3 kg)  
8-cell battery pack  
6 lbs (2.7 kg)  
16-cell battery pack

### SOUND LEVEL

39.9 dB(A) @ Setting 2

### BATTERY DURATION

8-cell battery  
up to 8 hours  
16-cell battery  
up to 16 hours

### POWER

AC Power: 100-240  
VAC, 50-60Hz  
DC Power: 11-18 VDC  
(10 max amp)  
Battery: Lithium Ion

### WARRANTY

2 years  
(Includes sieve beds)  
1 year on battery and  
power supplies



Caire Inc.  
2200, Airport Industrial Dr. Ste.  
500 Ball Ground, GA 30107 USA  
<https://www.caireinc.com/product>