

IntelliSense®

Trusted Brand



Mesh Nebulizer

Specifications

Nebulization Rate	**0.25 - 0.9 ml/min	Battery Life	Approx. 3 hours
Particle Size	*MMAD Approx. 5 µm MMAD = Mass Median Aerodynamic Diameter	Power Source	2 "AA" alkaline batteries 1.5 V
Medication Container Capacity	10 ml	Dimensions	*** Approx. 45 (W) x 125 (H) x 70 (D) mm
Residual Volume	Approx. 0.5 ml	Weight	*** Approx. 110g (excluding batteries)
Operation Mode	Continuous operation	Operating Temperature/ Humidity/Air Pressure	+10 °C to +40 C, 30 to 85 % RH (non-condensing), 800 to 1060 hPa
Sound	Less than 25 dB	Storage and Transport Temperature/Humidity	-20 °C to +60 °C, 10 to 90 % RH (non-condensing)
		Product	OMRON Mesh Nebulizer
		Model	NE-U200

* Measured values reflect internal data obtained by NEXT GENERATION IMPACTOR (NGI), measured with 2.5 % NaF aqueous solution, according to EN 13544-1:2007 + A1 :2009 and ISO27427:2010.
 ** The Nebulization rate is measured with saline 0.9% solution at 23 °C at spray angle of 30 ° with 3 ml of medication.
 *** Main unit assembled with medication container and disposable mesh.

MANUFACTURED BY: OMRON DALIAN CO., LTD.
 NO.28, DONG BEI ER STREET, ECONOMIC AND TECHNICAL DEVELOPMENT ZONE,
 DALIAN 116600, CHINA.

MADE IN CHINA

TO REPORT PRODUCT COMPLAINT DIAL TOLL FREE NO.: 1800-209-2505 OR EMAIL AT DSRM@LUPIN.COM

IMPORTED BY:
 OMRON HEALTHCARE INDIA PVT. LTD,
 6TH FLOOR, B-BLOCK, SEWA TOWER, PLOT NO-19,
 SECTOR-18, UDYOG VIHAR, MARUTI INDUSTRIAL COMPLEX , GURGAON, HARYANA (INDIA) - 122008

IMPORT LIC NO. IMP/MD/2021/000021

® REGD. TRADE MARK OF LUPIN LTD.



MARKETED BY :
LUPIN LTD.
 159, C.S.T. ROAD, KALINA,
 SANTACRUZ (E), MUMBAI - 400 098, INDIA.

Lot : **20210322U**
 Mfg. Date : **APR - 21**
 Imported on : **DEC - 21**
 Packed on : **NOV - 21**

MAXIMUM RETAIL PRICE : ₹8450/-
 (INCLUSIVE OF ALL TAXES)





LUPIN Nebster®
Mesh Nebulizer



FREE 30 PCS
LUPIN NEBSTER DISPOSABLE MESH



Disposable mesh (30 pcs)

LUPIN Nebster®
Disposable Mesh



Replaceable Mesh is available at your nearby Pharmacy/Chemist

Free sample with Lupin Nebster Mesh Nebulizer
For Lot, Mfg. Date & Imported on see carton.

RO

30 pcs

IntelliSense®
MOST RECOMMENDED BRAND BY DOCTORS

Most Trusted Brand

LUPIN
Nebster[®]

Mesh Nebulizer

EFFICIENT




EASY-TO-USE
DISPOSABLE MESH



PORTABLE



LOW NOISE



LUPIN
Nebster[®]
Mesh Nebulizer

Box Contents	Qty.
Carry Case	1 unit.
Main unit assembled with medication container	1 nos.
Adult Mask	1 nos.
Child Mask	1 nos.
Mouthpiece	1 nos.
2 AA Alkaline batteries	1 set.
Instruction Manual	1 nos.
Quick set-up Guide	1 nos.
Lupin Nebster Disposable Mesh*	30 pcs.

*Replaceable Mesh is available at your nearby Pharmacy/Chemist

User Instructions:

1. Lupin Nebster is a medical device, follow your Physician's instructions for dose and regimen of medication.
2. Device is not designed to function as an inhaler for moistening the mucous membranes of the nose and throat.

Medical Purpose – This product is intended to be used for inhaling medication for respiratory disorders.
Do not operate the nebulizer with the medication container empty.
Do not use the product for any purpose other than inhalation.
Do not block the air opening in the medication container during nebulization.
This device does not contain medicinal or biological substance.
Non sterile