

Oxymed i Series for OSA & COPD



TECHNICAL SPECIFICATIONS

Model Name	CPAP	BPAP
Model No	i-series C5	i-series P1
Modes	CPAP AUTO CPAP	CPAP, S, T, S/T, VAT, Pc
Pressure range	4-20cm H ₂ O	4-30cm H ₂ O
Pressure accuracy	± 0.2 cm H ₂ O	± 0.2 cm H ₂ O
Max. operation Pressure	30 cm H ₂ O	30 cm H ₂ O
COMF Pressure relief	1-3 levels	1-3 levels
Ramp time	0 to 45 min (5-minute increments)	0 to 45 min (5-minute increments)
Humidification Level	1-5 levels (113 to 185°F/23 to 85°C)	1-5 levels (113 to 185°F/23 to 85°C)
Rise Time	1-4 levels	1-6 levels (S,T,S/T)
Split night	YES	-
Data Storage Capacity	8G USB disk	8G USB disk
Weight	1.72 kg	1.72 kg
Mean Sound Level	≤ 30dB	≤ 30dB



Separable design




Preheat



M+ cloud platform




2 years Warranty*




Safe

AST

Automatic synchronization technology




Low-noise blower



VAT

Volume Assured Technology



Comfortable

Preheat

5 minutes from 23°C to 45 °C




One button control



COMF

Pressure relief at the end of exhalation



User-friendly

Separable design



USB storage



* Terms & Conditions apply