



Durable Nebulizer for Fast & Effective Relief

Reliable Respiratory Therapy with Medikart



Compressor Nebulizer for Hospital Use

A Reliable Method to Relieve Respiratory Symptoms



1 Year
Replacement
Warranty



Medikart Compressor Nebulizer MK NE-C05

Technical Specifications



Nebulization Rate:
>0.35 ml/min



Particle Size:
MMAD 2.5 μ m



Noise Level:
 \leq 55 dB(A)



Operating Pressure:
1.1 – 1.5 bar



Power Supply:
220V / 50Hz



Overheat Protection



What's in the Box



Medikart Compressor Nebulizer
MK NE-C05 Main Unit



Nebulizer Medication Cup



Adult Mask




Child Mask



Mouthpiece



Air Tubing Pipe



Air Filters



Power Cord



User Manual /
Instruction Guide



For Any Issue Contact Medikart



+91 76690 11322



+91 98119 57360



info@medikart.co.in



www.medikart.co.in



MediKart Healthcare Systems Pvt. Ltd.,
First Floor-45 Omaxe Square Jasola
District Centre,
Behind Apollo Hospital,
New Delhi – 110025.

Follow Us



Facebook



Instagram



Threads



YouTube



LinkedIn



Medikart Compressor Nebulizer MK NE-C05 – Hospital Use

Reliable Respiratory Therapy for Hospitals & Clinics



Product Overview

The Medikart Compressor Nebulizer MK NE-C05 is designed for hospital and clinical respiratory therapy, delivering consistent and efficient medication delivery for patients with breathing disorders.

Key Benefits



-  • Powerful compressor technology
-  • Fast and efficient nebulization
-  • Suitable for hospital and clinical use
-  • Continuous operation capability
-  • Easy to operate and maintain
-  • Durable and reliable design



Technical Specifications

Model	MK NE-C05
Type	Compressor Nebulizer
Application	Hospital / Clinical Use
Nebulization Rate	>0.35 ml/min
Particle Size	MMAD 2.5 µm
Noise Level	≤55 dB(A)
Operating Pressure	1.1 – 1.5 bar
Power Supply	220V/50Hz
Weight	approx 1.6 kg
Accessories Included	Mask (Adult & Child), Mouthpiece, Tubing, Medication Cup

Applications



• Asthma Treatment



• COPD Management



• Respiratory Therapy



• Hospital Patient Care