

**AGEasy**  
Max Group

# MEDICAL OXYGEN CONCENTRATOR

---

5L



## Description:

The AGEasy 5L Oxygen Concentrator is a reliable, modern solution designed to provide consistent supplemental oxygen with ease and confidence. Compact yet powerful, it delivers up to 5L per minute of oxygen flow, ensuring dependable performance for both home and healthcare use.

Built with a digital display, touch controls, and an oxygen purity indicator, the AGEasy concentrator is simple to operate, even for seniors. Its low noise level and lightweight, one-hand carry design make it ideal for daily use, blending advanced technology with user-friendly convenience.

Additionally, the concentrator combines efficiency, safety, and comfort to support better breathing and peace of mind. The in-built nebulizer function further enhances its utility, allowing seamless respiratory therapy without the need for a separate device.

## Key Features:

5 L/min Flow Range	Consistent oxygen delivery (0.5–5 L/min)
93% ±3% Oxygen Purity	High concentration for effective therapy
Digital Display & Touch Buttons	Modern, easy-to-use interface
Oxygen Purity Indicator	Real-time monitoring of oxygen quality
Room Temperature Display	Shows ambient temperature on-screen
Quiet Operation	Lower noise level for comfortable use
Compact Design	Lightweight design, easy to carry with one hand
Durable Build	Made for everyday use at home or in healthcare settings
In-built Nebulizer	Dual functionality for oxygen therapy and medication delivery



## Patient's Advantages:

- Energy efficient: only 280 W power consumption
- Auto shut off and safety protections
- Clear, bright display for seniors and caregivers
- Easy maintenance with accessible components

## Specifications:

Oxygen Flow	0.5 – 5 L/min
Power	280 W
Voltage/Frequency	220 V, 50 Hz
Dimensions	35 × 25.4 × 55.4 cm
Gross Weight	13.9 kg
Warranty	2 years / 10,000 hours

## What's in the Box?

- AGEasy 5L Oxygen Concentrator
- Nebulizer Kit
- Humidifier Bottle
- Humidifier Bottle Connector
- Nasal Cannula
- Dust Filter (Extra)
- HEPA Filter (Extra)
- Power Cord
- User Manual







**Senior Care Starts With Us**


**AGEasy**  
 **Max Group**

**Contact:**

---

 +91-97173 05336

 [www.ageasy.co.in](http://www.ageasy.co.in)

 2nd Floor, Plot No. 65, Sector 44, Gurugram, Haryana 122003