

# F3 Fetal Monitor



# F3 Fetal Monitor



Rolling Stand



Wall Mount

Based on our own design language, F3 gets an iF design award for its compact lightweight and ergonomic design.


With its elegant design, bright color screen, and enhanced functionality, F3 fetal monitor would address the needs of obstetrical departments in both clinics and hospitals. F3 is amazingly lightweight and portable, which makes it suitable for outpatient service. With the latest technologies from EDAN, F3 offers an extensive set of monitoring parameters such as FHR, TOCO, DECG, IUP and fetal movement.

## Multi-display Modes



Switching among three display modes, user may choose the most suitable one for clinical use.

## Simple Workflow




Functional keys and knob provide various shortcuts to achieve functions for clinical use. The 'start' button could be configured to integrate patient info and printing. This could help the doctor to simplify the workflow and work with 1 button only.

## USB Storage Capacity



F3 supports USB interface, providing enlarged storage capacity, which is a quite good solution for house call service.

## Data Management



With the help of EDAN Insight software, doctors can transmit CTG data from F3 to PC for data management.

### Convenient to Carry or Transport

- Compact, portable and lightweight
- Long life Li-ion battery for outpatient service

### Easy to Read & Operate

- 5.6" foldable color TFT screen
- Large numeric and waveform display
- Simple workflow with silicone buttons for operations

### Advanced FHR Detection Technology

- Signals Overlap Verification to differentiate twins FHR
- FHR signal quality indicator helps optimize the probe position

### Flexible Printing Options

- Built-in thermo printer
- Print on 150/152 mm width recording paper
- 1/2/3 cm/min real-time printing speed
- 15 mm/sec fast printing speed for history traces

### Powerful Data Management

- 60 hours built-up memory for seamless monitoring
- MFM-CNS network for remote monitoring
- MFM-CNS Lite for data management on PC
- USB port to enlarge storage capacity



**F3**  
Fetal Monitor

**Global Headquarters:**

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen  
518122 P.R. China | +86.755.26898326 | [www.edan.com](http://www.edan.com) | [info@edan.com](mailto:info@edan.com)

**U.S. and Canada inquiries:**

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126  
+1.858.750.3066 | [www.edandiagnostics.com](http://www.edandiagnostics.com) | [edan-info@edandiagnostics.com](mailto:edan-info@edandiagnostics.com)

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.