



ALWAYS READY
TO SAVE LIVES

FLOOR STANDING AED CABINETS

WITH AUDIO VISUAL ALARM BEACON SYSTEM

High visibility. Instant access. Maximum safety. Our floor standing AED Cabinets are designed to protect life-saving devices and ensure quick accessibility in critical moments.

KEY FEATURES



AUDIO VISUAL ALARM BEACON

High intensity flashing light and loud alarm for maximum visibility and alert.



STRONG & DURABLE CONSTRUCTION

Made from heavy-duty CRCA steel with anti-corrosion powder coating for long life.



SECURE STORAGE

Lockable door with toughened glass window to protect AED from theft, tampering & dust.



HIGH VISIBILITY DESIGN

Bright orange color with universal AED signage for quick identification in emergencies.



INSTANT ACCESS

Ergonomically designed height ensures AED is easily accessible when every second counts.



EASY INSTALLATION

Floor standing design for stable placement in public spaces, offices, hospitals, schools, malls & industries.



Open View



Closed View

SPECIFICATIONS

Cabinet Type	Floor Standing
Material	CRCA Steel
Finish	Powder Coated (Bright Orange)
Door	Lockable with Toughened Glass Window
Alarm System	Audio Visual Alarm Beacon (110 dB)
Power Supply	220V AC / 12V DC (Adapter)
Compatibility	Universal for most AED Brands
Dimensions (H x W x D)	1350 x 400 x 300 mm (Approx.)
Applications	Indoor & Semi-Outdoor Use

IDEAL FOR



OFFICES



HOSPITALS



SCHOOLS &
COLLEGES



MALLS &
PUBLIC PLACES



INDUSTRIES &
WORKPLACES



PROTECTS LIFE
Saving Devices



QUICK ALERT
System



DURABLE &
RELIABLE



LOW
MAINTENANCE



MEETS SAFETY
STANDARDS