

OXYGEN CONCENTRATOR



UNO 10 Oxygen Concentrator

UNO Supply introduces the latest product in the Oxygen Concentrators category. UNO10. An innovative and evolved product. It is the result of our pursuit to provide the absolute best to our customers.

UNO 10 is lighter and more compact, with more functionality and benefits. It guarantees complete user satisfaction and convenience.

For medical use in hospitals, homes, and remote areas where oxygen supply is limited. Increases oxygen concentration for a single patient.

Contraindication: Patients with severe carbon monoxide poisoning.

A portable oxygen concentrator using Pressure Swing Adsorption (PSA) with zeolite molecular sieve to separate oxygen from air. Produces 90–96% oxygen continuously at normal temperature and pressure. Does not affect surrounding air oxygen levels.



SPECIFICATIONS

Delivery Rate (LPM)	0 – 10 (LPM)
Maximum Flow (LPM)	10 (LPM)
Purity (%)	93 ± 3 (2–10 L)
Power Consumption (Watts)	600
Weight (Kg)	24 Kgs
Dimensions (LxWxH in mm)	372 × 340 × 612mm
OCI	Yes
Sound Level	≤ 50dB
Alarms	Power failure, Overheat, Oxygen low purity
Dual Flow	Optional
Electrical Requirement	230V ; 50Hz



Head Office: 238B, A.J.C. Bose Road,
Kolkata 700 020
E-mail: info@unosupply.com
Customer Care: +91 98300 56776
Website: www.unosupply.com