

# BREATHWEAR MASK SERIES

## YF-02



Fit



Comfort



Convenience



Light & Soft

- Type: Full Face Mask
- Therapy Pressure Range :4-30 cmH2O
- Tube Connection :  $\varphi 22\text{mm}$
- Flow Resistance(At 50 L/min):  $\leq 0.5\text{cmH}_2\text{O}$
- Flow Resistance(At 100 L/min):  $\leq 1.5\text{cmH}_2\text{O}$
- Expiratory Resistance with Anti Asphyxia Valve Open (At 50 L/min):  $\leq 2.0\text{cmH}_2\text{O}/(\text{L/S})$
- Inspiratory Resistance with Anti Asphyxia Valve Open (At 50 L/min):  $\leq 2.0\text{cmH}_2\text{O}/(\text{L/S})$
- Anti Asphyxia Valve Open Pressure:  $\leq 2.0\text{cmH}_2\text{O}/(\text{L/S})$
- Anti Asphyxia Valve Closed Pressure:  $\leq 2.0\text{cmH}_2\text{O}/(\text{L/S})$



## PACKAGING PARAMETERS

Model	QTN (PCS/CTN)	G.W/N.W (KG)	Size (cm)
YF-02	10	2/1.6	35×28×28

**JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD**

Sales Lenter Address: Huanyuan tast Road No. I, Xuzhuang Sortware Park, Nanjing

Jianasu Province, PR. China

Tel: 86-25-8713 6530 Fax: 86-25-8713 6530

E-mail: yuwellsales@yuyue.com.cn ; sales@yuyue.com.cn

www.yuwell.com